

# Spot Stain Removal

# RS

## Red Spot



212



PRODUCT CODE	212
APPEARANCE	clear to slightly hazy pale straw liquid
REACTION	slightly alkaline
SPECIFIC GRAVITY	1.020 – 1.040 at 20°C
STORAGE	above 5°C
STANDARD PACK SIZE	5L
ALSO AVAILABLE IN	10ml roll-on applicator, 250ml, 20L

## Professional stain spotter for treating blood, milk, food & body fluid stains

Developed in consultation with professional cleaners for tackling specific problem stains.

Treat blood, milk, cream, egg, mayonnaise and chocolate stains. Remove urine, vomit and body fluids.

Ready to use straight from the bottle.

Available as 10ml mini roll-on or 250ml applicator bottles. Refills available in 5L drums.

### How to use

Refer to Ideal Spotter Chart.

Apply spotter directly to stain. Allow time to penetrate. Gently work stain with a spotting bone, brush or Ideal roll-on applicator. Flush area with water or steam. Avoid dry steam on blood stains. For set blood stains saturate with red spot and allow extended time for solubilisation. Clean treated garments as normal.

Dyestuffs can often display poor fastness to spotting, so we recommend that seam tests are always carried out before application.

It's always better to clean as soon as possible after application.

You should only pre spot garments in time for the next load.

For further process information speak to an Ideal technical representative.

### Health & safety

Available on request. Store above 5°C.



Ideal Manufacturing Ltd.

Atlas House, Burton Road, Finedon, Northants, NN9 5HX

T 01933 681616

hello@idealmanufacturing.com

Follow us @IdealManLtd

www.idealmanufacturing.com

# BRIGHTER LAUNDRY

